



Release Your *Full Potential*

Neoss® Guide Kit – Flexible



Intelligent Simplicity





One guided surgery kit for all three Neoss systems Newly designed drilling protocol

All our implant systems can be inserted with a single set of drills and a single guide kit, reducing the complexity of surgical preparation and simplifying your inventory.

Full depth flexibility for specific clinical need

Streamlined workflow

The Flexible Kit setup is designed to work with the different offsets of G3, G5 and G7 (3mm, 5mm and 7mm). This allows for a guided-surgery design that suits different clinical needs. By utilising the three-line markings on the tools and combinations of hubs and implant length, an experienced clinician can control drilling depths with the same set of tools.




The new sleeve design enhances seating stability and supports precise orientation control. Open sleeves reduce the vertical height requirement, improving access in challenging cases.

	Closed Sleeve	Open Sleeve	Implant Diameter
Narrow – N (Ø4.2)			Ø3.25 – Ø3.5
Wide – W (Ø5.5)			Ø3.5 – Ø5.0

Achieve precise *pre-designed* implant orientation Full digital workflow compatibility

Integrated orientation control on the implant inserter and sleeves, which is also fully reflected in our digital design databases, enables accurate compatibility with pre-designed and prefabricated temporary restorations for immediate loading cases.

Design databases are available in leading planning software programs, with more on the way.



Assortment

Neoss® Guide Kit Ø3.25 – Ø5.0 Flexible

Art. No.	Description
----------	-------------

51258	Neoss® Guide Kit Ø3.25 – Ø5.0 G3/5/7 Flexible
51260	Neoss® Guide Instrument Kit (refer to next page)


51258 includes:

Art. No.	Description
----------	-------------

51201	Neoss® Guide - Punch - Narrow
51202	Neoss® Guide - Mill - Narrow G3/5/7
51203	Neoss® Guide - Start Drill - Narrow
51204	Neoss® Guide - Punch - Wide
51205	Neoss® Guide - Mill - Wide G3/5/7
51206	Neoss® Guide - Start Drill - Wide
51207	Neoss® Guide - Drill Hub I, Narrow
51208	Neoss® Guide - Drill Hub II, Narrow
51209	Neoss® Guide - Drill Hub III, Narrow
51210	Neoss® Guide - Drill Hub IV, Narrow
51211	Neoss® Guide - Drill Hub V, Narrow
51212	Neoss® Guide - Drill Hub I, Wide
51213	Neoss® Guide - Drill Hub II, Wide
51214	Neoss® Guide - Drill Hub III, Wide
51215	Neoss® Guide - Drill Hub IV, Wide
51216	Neoss® Guide - Drill Hub V, Wide
51217	Neoss® Guide - Twist Drill, Ø2.2 × 7-15 mm
51218	Neoss® Guide - Twist Drill, Ø2.85 × 7-15 mm
51219	Neoss® Guide - Twist Drill, Ø3.0 × 7-15 mm
51220	Neoss® Guide - Twist Drill, Ø3.2 × 7-15 mm
51221	Neoss® Guide - Twist Drill, Ø3.4 × 7-15 mm
51222	Neoss® Guide - Twist Drill, Ø3.6 × 7-15 mm
51223	Neoss® Guide - Twist Drill, Ø3.9 × 7-15 mm
51224	Neoss® Guide - Twist Drill, Ø4.1 × 7-15 mm
51225	Neoss® Guide - Twist Drill, Ø4.4 × 7-15 mm
51226	Neoss® Guide - Twist Drill, Ø4.6 × 7-15 mm
51227	Neoss® Guide - Countersink, Ø3.25 S - Narrow G3/5/7
51228	Neoss® Guide - Countersink, Ø3.5 S/T/E - Narrow G3/5/7
51248	Neoss® Guide - Countersink, Ø3.5 S/T/E - Wide G3/5/7
51229	Neoss® Guide - Countersink, Ø4.0 S - Wide G3/5/7
51230	Neoss® Guide - Countersink, Ø4.5 S - Wide G3/5/7
51231	Neoss® Guide - Countersink, Ø5.0 S - Wide G3/5/7
51232	Neoss® Guide - Countersink, Ø4.0 T/E - Wide G3/5/7
51233	Neoss® Guide - Countersink, Ø4.5 T/E - Wide G3/5/7
51234	Neoss® Guide - Countersink, Ø5.0 T/E - Wide G3/5/7
51238	Neoss® Guide - Implant Inserter, SP, Wide
51239	Neoss® Guide - Implant Inserter, SP, Narrow
51240	Neoss® Guide - Implant Inserter, NP, Narrow
51241	Neoss® Guide Kit System Box




Neoss Guide Kit



51260 includes:

Art. No.	Description
51121	Ratchet
51144	Wrench Adapter
51137	Implant inserter, SP 17mm
51126	Manual Handle
51139	Neo Screwdriver, Machine 22mm



Guide Sleeves

Guide Sleeves can be ordered separately for Guide.

Art. No.	Description
51242	Neoss® Guide - Open Sleeve, Narrow - 5 pcs
51243	Neoss® Guide - Closed Sleeve, Narrow - 5 pcs
51244	Neoss® Guide - Open Sleeve, Wide - 5 pcs
51245	Neoss® Guide - Closed Sleeve, Wide - 5 pcs



Implant Mounts (Optional)

Implant Mounts are optional and can be ordered separately.

Art. No.	Description
51235	Neoss® Guide - Implant Mount, NP, Narrow G5
51236	Neoss® Guide - Implant Mount, SP, Narrow G5
51237	Neoss® Guide - Implant Mount, SP, Wide G5



You can purchase 51259 Neoss® Guide add-on-kit Countersink and Mill G5 for extra components in hard stop mode. Please contract your local Neoss Representatives for the options.

AUSTRALIA & NEW ZEALAND
Neoss Australia Pty. Ltd.
 G04 Ground Floor North Tower, John Oxley Centre
 339 Coronation Drive
 Milton Queensland 4064
 Australia
 T: + +61 7 3216 0165
 E: info.au@neoss.com

Order your products
 24/7 in our Online Stores



Australia
 at shop-au.neoss.com



New Zealand
 at shop-nz.neoss.com

For our Terms & Conditions consult your local website on www.neoss.com

Contact your local representative for product availability details.
 All products are not available in every market.

17519_0 EN 26-03 © Neoss Group, 2026. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss Group. Neoss and the N logo are registered trademarks of Neoss Group. NeoScan™, NeoPro™, and NeoConnect™ are trademarks of Neoss Group.