



25 Years
Anniversary

Product Catalogue

2026-27

Intelligent Simplicity

Explore Intelligent Simplicity

Neoss offers intelligent solutions that are intuitively easy to use. Our products allow dental professionals to provide reliable and cost-effective treatments to their patients with predictable long-term results. With ingenuity and ambition, we strive to set new standards, leading the market with integrity. In developing smart treatment solutions and working closely with our customers, Neoss makes the complex less complicated. We call that Intelligent Simplicity.



Implants

- 4 **Neoss ProActive® Edge Implants**
- 6 **Neoss ProActive® Tapered Implants**
Neoss ProActive® Wide Implants
Neoss ProActive® Sinus Implants
- 7 **Neoss ProActive® Straight Implants**
- 8 **Neoss ProActive® Straight NP Ø3.25 Implants**
Cover Screw
- 9 **Healing Abutments**
Esthetic Healing Abutments with ScanPeg
Healing Abutments
- 10 **Provisional Abutments**
Esthetic Tissue Formers
Provisional Titanium Abutments

Instruments & More

- 11 **Drills**
- 13 **Kits**
Surgical & Prosthetic Trays
Surgical & Prosthetic Kits with Instruments and Drills
- 14 **Neoss® Guide Kits – Hard Stop & Flexible** **NEW**
- 15 **Kit-Specific Components & Assortments**
- 16 **Instruments**
Demo Kit

Digital

- 17 **Resonance Frequency Device (RFA)**
NeoTell®
Intraoral Scanners
NeoScan™ 1000
NeoScan™ 2000
- 18 **Digital Impression**
Esthetic Healing Abutments with ScanPeg
ScanPeg
Scan Bodies
ScanKeys
Digital Analogs
- 19 **Digital Prosthetics**
NeoBase® Abutments SSC (Straight Screw Channel)*
NeoBase® Abutments ASC (Angulated Screw Channel)*
Neoss TiBase Abutments (compatible with Dentsply Sirona®)
Ti NeoLink®

Prosthetics

- 20 **Impressions & Replicas**
Impression Copings
Replicas
- 21 **Final Abutments**
Prepable Titanium Abutments
Prepable Titanium Abutments, Alternative Emergence Profiles
Zirconia Abutments
- 22 **Final Restorations – Casting**
Gold NeoLink® (Direct Casting)
CoCr Abutments (Indirect Casting)
Burnout Abutments (Direct Casting)
Neoss® Medentika Blanks **NEW**
- 23 **Neoss4+ treatment solution**
Multi-Unit Abutments
- 24 **Access Abutments**
Multi-Unit/Access Screws
- 25 **Overdenture Abutments**
Equator Abutments
Screw Guide
Screws
Overview Abutment, Prosthetic and Laboratory Screws

26–27 ARC Solutions

Biomaterials

- 28 **Bone Graft Materials**
NeoGenix XP
PTFE Sutures
Collagen Membranes
- 29 **Non-Resorbable Membranes**
Ti-Reinforced PTFE Membranes
NeoGen® Cape
Non-Reinforced PTFE Membranes
Fixation Kit

Drilling and Drill Depth Guides

- 30 **ProActive Edge Implants**
ProActive Tapered Implants
- 31 **ProActive Straight Implants**
Drill Depth Guides
- 32 **ProActive Wide and Sinus Implants**

Clinical Organisers – Surgical Set-Up



- 32 **ProActive Edge Implants**
- 33 **ProActive Tapered Implants**
ProActive Straight Implants

Support

- 34 **Prosthetic Solution Standard and Narrow Platform, overview**
Torque and Speed recommendation, overview
- 35 **Symbols, names and abbreviations**
You can now buy your products online!
Disclaimer of Liability

Neoss ProActive® implants – *Proven predictability* in implant treatment since 2002

All implants produced of highest quality, Commercially Pure Titanium Grade 4.



Neoss Implant

The Implant includes Implant and Cover Screw.





Neoss ProActive® Edge Implants

Art. No.	Description	A	B
21417	Implant – Neoss ProActive® Edge Ø3.5×9 mm	4.0mm	1.6mm
21418	Implant – Neoss ProActive® Edge Ø3.5×11mm		
21419	Implant – Neoss ProActive® Edge Ø3.5×13mm		
21420	Implant – Neoss ProActive® Edge Ø3.5×15mm		
21421	Implant – Neoss ProActive® Edge Ø3.5×17mm		
21423	Implant – Neoss ProActive® Edge Ø4.0×9 mm	4.3mm	1.9mm
21424	Implant – Neoss ProActive® Edge Ø4.0×11mm		
21425	Implant – Neoss ProActive® Edge Ø4.0×13 mm		
21426	Implant – Neoss ProActive® Edge Ø4.0×15 mm		
21427	Implant – Neoss ProActive® Edge Ø4.0×17 mm		
21429	Implant – Neoss ProActive® Edge Ø4.5×9 mm	4.9mm	2.2mm
21430	Implant – Neoss ProActive® Edge Ø4.5×11mm		
21431	Implant – Neoss ProActive® Edge Ø4.5×13 mm		
21432	Implant – Neoss ProActive® Edge Ø4.5×15 mm		
21433	Implant – Neoss ProActive® Edge Ø4.5×17 mm		
21435	Implant – Neoss ProActive® Edge Ø5.0×9 mm	5.4mm	2.4mm
21436	Implant – Neoss ProActive® Edge Ø5.0×11mm		
21437	Implant – Neoss ProActive® Edge Ø5.0×13 mm		
21438	Implant – Neoss ProActive® Edge Ø5.0×15 mm		
21439	Implant – Neoss ProActive® Edge Ø5.0×17 mm		



SP = Standard Platform
Ø3.5–6.5 Implants



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments in PEEK in two different sizes and Healing Abutment Screw.

Neoss ProActive® Edge Implants

Art. No.	Description	A	B
21301	Implant Kit – Neoss ProActive® Edge Ø3.5×9 mm		
21302	Implant Kit – Neoss ProActive® Edge Ø3.5×11 mm	4.0 mm	1.6 mm
21303	Implant Kit – Neoss ProActive® Edge Ø3.5×13 mm		
21304	Implant Kit – Neoss ProActive® Edge Ø3.5×15 mm		
21305	Implant Kit – Neoss ProActive® Edge Ø4.0×9 mm		
21306	Implant Kit – Neoss ProActive® Edge Ø4.0×11 mm	4.3 mm	1.9 mm
21307	Implant Kit – Neoss ProActive® Edge Ø4.0×13 mm		
21308	Implant Kit – Neoss ProActive® Edge Ø4.0×15 mm		
21309	Implant Kit – Neoss ProActive® Edge Ø4.5×9 mm		
21310	Implant Kit – Neoss ProActive® Edge Ø4.5×11 mm	4.9 mm	2.2 mm
21311	Implant Kit – Neoss ProActive® Edge Ø4.5×13 mm		
21312	Implant Kit – Neoss ProActive® Edge Ø4.5×15 mm		
21313	Implant Kit – Neoss ProActive® Edge Ø5.0×9 mm		
21314	Implant Kit – Neoss ProActive® Edge Ø5.0×11 mm	5.4 mm	2.4 mm
21315	Implant Kit – Neoss ProActive® Edge Ø5.0×13 mm		
21316	Implant Kit – Neoss ProActive® Edge Ø5.0×15 mm		



ProActive Edge Ø3.5




ProActive Edge Ø4.0



ProActive Edge Ø4.5



ProActive Edge Ø5.0



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0 – 5.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.



Neoss ProActive® Tapered Implants

Art. No.	Description	A	B
21221	Implant Kit – Neoss ProActive® Tapered Ø3.5×9mm		
21222	Implant Kit – Neoss ProActive® Tapered Ø3.5×11mm	4.0mm	2.1mm
21223	Implant Kit – Neoss ProActive® Tapered Ø3.5×13mm		
21224	Implant Kit – Neoss ProActive® Tapered Ø3.5×15mm		
21227	Implant Kit – Neoss ProActive® Tapered Ø4.0×9mm		
21228	Implant Kit – Neoss ProActive® Tapered Ø4.0×11mm	4.3mm	2.3mm
21229	Implant Kit – Neoss ProActive® Tapered Ø4.0×13mm		
21230	Implant Kit – Neoss ProActive® Tapered Ø4.0×15mm		
21233	Implant Kit – Neoss ProActive® Tapered Ø4.5×9mm		
21234	Implant Kit – Neoss ProActive® Tapered Ø4.5×11mm	4.9mm	2.8mm
21235	Implant Kit – Neoss ProActive® Tapered Ø4.5×13mm		
21236	Implant Kit – Neoss ProActive® Tapered Ø4.5×15mm		
21239	Implant Kit – Neoss ProActive® Tapered Ø5.0×9mm		
21240	Implant Kit – Neoss ProActive® Tapered Ø5.0×11mm	5.4mm	3.3mm
21241	Implant Kit – Neoss ProActive® Tapered Ø5.0×13mm		
21242	Implant Kit – Neoss ProActive® Tapered Ø5.0×15mm		
21245	Implant Kit – Neoss ProActive® Tapered Ø5.5×9mm		
21246	Implant Kit – Neoss ProActive® Tapered Ø5.5×11mm	5.9mm	3.8mm
21247	Implant Kit – Neoss ProActive® Tapered Ø5.5×13mm		

Neoss ProActive® Wide Implants



Art. No.	Description	A	B
21252	Implant Kit – Neoss ProActive® Ø6.0×7mm		4.5mm
21250	Implant Kit – Neoss ProActive® Ø6.0×9mm	6.0mm	4.3mm
21251	Implant Kit – Neoss ProActive® Ø6.0×11mm		4.3mm

Neoss ProActive® Sinus Implants

Art. No.	Description	A
21211	Implant Kit – Neoss ProActive® Ø6.5×7mm	
21212	Implant Kit – Neoss ProActive® Ø6.5×9mm	6.7mm
21213	Implant Kit – Neoss ProActive® Ø6.5×11mm	



SP = Standard Platform
Ø3.5 – 6.5 Implants



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0 – 6.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.

Neoss ProActive® Straight Implants

Art. No.	Description	A	B
21181	Implant Kit – Neoss ProActive® Straight Ø3.5×7 mm		
21182	Implant Kit – Neoss ProActive® Straight Ø3.5×9 mm		
21183	Implant Kit – Neoss ProActive® Straight Ø3.5×11 mm	4.0 mm	2.4 mm
21184	Implant Kit – Neoss ProActive® Straight Ø3.5×13 mm		
21185	Implant Kit – Neoss ProActive® Straight Ø3.5×15 mm		
21186	Implant Kit – Neoss ProActive® Straight Ø3.5×17 mm		
21187	Implant Kit – Neoss ProActive® Straight Ø4.0×7 mm		
21188	Implant Kit – Neoss ProActive® Straight Ø4.0×9 mm	4.0 mm	2.8 mm
21189	Implant Kit – Neoss ProActive® Straight Ø4.0×11 mm		
21190	Implant Kit – Neoss ProActive® Straight Ø4.0×13 mm		
21191	Implant Kit – Neoss ProActive® Straight Ø4.0×15 mm		
21192	Implant Kit – Neoss ProActive® Straight Ø4.0×17 mm		
21193	Implant Kit – Neoss ProActive® Straight Ø4.5×7 mm		
21194	Implant Kit – Neoss ProActive® Straight Ø4.5×9 mm	4.5 mm	3.2 mm
21195	Implant Kit – Neoss ProActive® Straight Ø4.5×11 mm		
21196	Implant Kit – Neoss ProActive® Straight Ø4.5×13 mm		
21197	Implant Kit – Neoss ProActive® Straight Ø4.5×15 mm		
21198	Implant Kit – Neoss ProActive® Straight Ø4.5×17 mm		
21199	Implant Kit – Neoss ProActive® Straight Ø5.0×7 mm		
21200	Implant Kit – Neoss ProActive® Straight Ø5.0×9 mm	5.0 mm	3.7 mm
21201	Implant Kit – Neoss ProActive® Straight Ø5.0×11 mm		
21202	Implant Kit – Neoss ProActive® Straight Ø5.0×13 mm		
21203	Implant Kit – Neoss ProActive® Straight Ø5.0×15 mm		
21205	Implant Kit – Neoss ProActive® Straight Ø5.5×7 mm		
21206	Implant Kit – Neoss ProActive® Straight Ø5.5×9 mm	5.5 mm	4.2 mm
21207	Implant Kit – Neoss ProActive® Straight Ø5.5×11 mm		
21208	Implant Kit – Neoss ProActive® Straight Ø5.5×13 mm		



ProActive Straight Ø3.5



ProActive Straight Ø4.0



ProActive Straight Ø4.5



ProActive Straight Ø5.0




ProActive Straight Ø5.5



Neoss Implant Kit

The Neoss ProActive® NP Ø3.25 implant kit includes Implant, Cover Screw, Healing Abutment and Healing Abutment Screw.



ProActive Straight NP Ø3.25

Neoss ProActive® Straight NP Ø3.25 Implants


Art. No. Description		A	B
21176	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×9 mm	3.5mm	2.2mm
21177	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×11 mm		
21178	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×13 mm		
21179	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×15 mm		



Cover Screw

Cover Screw

Art. No. Description	
31408	Cover Screw SP
31409	Cover Screw NP



Healing Abutments

Esthetic Healing Abutments with ScanPeg


All Esthetic Healing Abutments include Titanium Screw and a ScanPeg


Art. No.	Description	A	B	C
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor	5mm	5.7mm	6.1mm
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor		5.2mm	5mm
31362	Esthetic Healing Abutment with ScanPeg – Canine		6.3mm	5.4mm
31363	Esthetic Healing Abutment with ScanPeg – Pre-molar		6.8mm	5.8mm
31364	Esthetic Healing Abutment with ScanPeg – Molar		7.8mm	7.3mm



Healing Abutments

Art. No.	Description	A	B
31159	Healing Abutment Ti 2mm	2mm	4mm
31160	Healing Abutment Ti 4mm	4mm	
31161	Healing Abutment Ti 6mm	6mm	
31162	Healing Abutment Ti 8mm	8mm	
31163	Healing Abutment Ti 10mm	10mm	
31148	Healing Abutment PEEK 5mm with Screw	5mm	5mm
31320	Healing Abutment PEEK 5mm with Screw, NP		4mm
31321	Wide Healing Abutment PEEK 3.3mm with Screw	3.3mm	6.5mm
31335	Wide Healing Abutment PEEK 5mm with Screw	5mm	6.5mm
31410	Healing Abutment PEEK 2.7 mm with Screw	2.7mm	5mm






Provisional Abutments

Esthetic Tissue Formers

All Esthetic Tissue Formers include Provisional Screw



Esthetic Tissue Formers

Art. No.	Description	A	B	C	D
31340	Esthetic Tissue Former – Wide Incisor	10 mm	4 mm	5.7 mm	6.1 mm
31341	Esthetic Tissue Former – Narrow Incisor			5.2 mm	5 mm
31342	Esthetic Tissue Former – Canine			6.3 mm	5.3 mm
31343	Esthetic Tissue Former – Pre-molar			6.8 mm	5.8 mm
31344	Esthetic Tissue Former – Molar			7.8 mm	7.3 mm



Provisional Ti Abutment Mono



Provisional Ti Abutment Multi



Provisional Ti Abutment Mono – Cement Retained



Provisional Ti Abutment, NP



Coping

20 mm

Provisional Titanium Abutments

All Provisional Abutments include Provisional Screw and Coping

Art. No.	Description	A	B	C
31267	Provisional Ti Abutment Mono	15.7 mm	0.7 mm	4.5 mm
31268	Provisional Ti Abutment Multi			
31332	Provisional Ti Abutment Mono – Cement Retained	10.7 mm	0.7 mm	3.8 mm
31241	Provisional Ti Abutment, NP			

Includes Provisional Titanium Abutment, Provisional Screw NP and Neo Laboratory Screw NP

The instruments *you need*

Drills

Art. No. Description

41153	Round Bur Ø1.8mm	
41152	Lance Drill Ø1.5×3/5/7mm	
41170	Pilot Drill Ø2.2/3.0mm	
41154	Twist Drill, Straight Ø2.2×7 – 17mm	
41156	Twist Drill, Straight Ø2.85×7 – 17mm	
41157	Twist Drill, Straight Ø3.0×7 – 17mm	
41158	Twist Drill, Straight Ø3.2×7 – 17mm	
41159	Twist Drill, Straight Ø3.4×7 – 17mm	
41160	Twist Drill, Straight Ø3.6×7 – 17mm	
41161	Twist Drill, Straight Ø3.9×7 – 17mm	
41162	Twist Drill, Straight Ø4.1×7 – 17mm	
41163	Twist Drill, Straight Ø4.4×7 – 17mm	
41164	Twist Drill, Straight Ø4.6×7 – 17mm	
41165	Twist Drill, Straight Ø4.9×7 – 17mm	
41166	Twist Drill, Straight Ø5.1×7 – 17mm	
41205	Twist Drill, Straight Ø6.0×7 – 11 (13)mm	
41193	Drill Kit, Straight Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9 and 4.4</i>	
41172	Twist Drill, Tapered Ø3.0×7 – 17mm	1.9mm
41173	Twist Drill, Tapered Ø3.4×7 – 17mm	2.3mm
41174	Twist Drill, Tapered Ø3.9×7 – 17mm	2.8mm
41175	Twist Drill, Tapered Ø4.4×7 – 17mm	3.3mm
41176	Twist Drill, Tapered Ø4.9×7 – 17mm	3.8mm
41206	Twist Drill, Tapered Ø5.5×7 – 17mm	4.3mm
41204	Twist Drill, Tapered Ø5.5×7 – 11 (13)mm	4.3mm
41192	Drill Kit, Tapered Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2 S, 3.0 T, 3.4 T, 3.9 T and 4.4 T</i>	
41177	Drill Kit, Edge Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2 S, 3.0 T, 3.4 T, 3.9 T and 4.4 T</i>	
41142	Drill Stop 7mm Ø2.2, 3.0, 3.4, 3.9 and 4.4mm	
41143	Drill Stop 9mm Ø2.2, 3.0, 3.4, 3.9 and 4.4mm	
41144	Drill Stop 11mm Ø2.2, 3.0, 3.4, 3.9 and 4.4mm	
41145	Drill Stop 13mm Ø2.2, 3.0, 3.4, 3.9 and 4.4mm	
41146	Drill Stop 15mm Ø2.2, 3.0, 3.4, 3.9 and 4.4mm	
41215	Twist Drill, Straight Ø2.2×7 – 17mm	
41216	Twist Drill, Straight Ø2.85×7 – 17mm	
41217	Twist Drill, Straight Ø3.0×7 – 17mm	
41218	Twist Drill, Straight Ø3.2×7 – 17mm	
41219	Twist Drill, Straight Ø3.4×7 – 17mm	
41220	Twist Drill, Straight Ø3.6×7 – 17mm	
41221	Twist Drill, Straight Ø3.9×7 – 17mm	
41222	Twist Drill, Straight Ø4.1×7 – 17mm	
41223	Twist Drill, Straight Ø4.4×7 – 17mm	
41224	Twist Drill, Straight Ø4.6×7 – 17mm	
41225	Twist Drill, Straight Ø4.9×7 – 17mm	
41226	Twist Drill, Straight Ø5.1×7 – 17mm	



Round Bur



Lance Drill



Pilot Drill



Twist Drill, Straight



Twist Drill, Tapered



Drill Stop

Single Use

Multiple Use



Twist Drill, Straight – Short



Twist Drill, Tapered – Short



Countersink



Countersink, Straight



Countersink, Tapered



Countersink, Edge



Bone Profiler



Bone Remover Cover Screw

Art. No. Description

41233 Drill Kit, Straight Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9 and 4.4

41184 Drill Kit, Straight – Short 7 – 13 mm
Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9, 4.4 and 4.9

41227 Twist Drill, Tapered Ø3.0×7 – 17 mm

41228 Twist Drill, Tapered Ø3.4×7 – 17 mm

41229 Twist Drill, Tapered Ø3.9×7 – 17 mm

41230 Twist Drill, Tapered Ø4.4×7 – 17 mm

41231 Twist Drill, Tapered Ø4.9×7 – 17 mm

41232 Twist Drill, Tapered Ø5.5×7 – 17 mm

41235 Drill Kit, Tapered Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T

41191 Drill Kit Tapered – Short 7 – 13 mm
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T, 4.4T, and 4.9T

41236 Drill Kit, Edge Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T

41139 Countersink, Straight Ø3.25

41110 Countersink, Straight Ø3.5

41111 Countersink, Straight Ø4.0

41112 Countersink, Straight Ø4.5

41133 Countersink, Straight Ø5.0

41127 Countersink, Straight Ø5.5

41178 Countersink, Tapered Ø3.5

41179 Countersink, Tapered Ø4.0

41180 Countersink, Tapered Ø4.5

41181 Countersink, Tapered Ø5.0

41182 Countersink, Tapered Ø5.5

41210 Countersink, Edge Ø3.5

41211 Countersink, Edge Ø4.0

41212 Countersink, Edge Ø4.5

41213 Countersink, Edge Ø5.0

41202 Countersink Ø6.0

41201 Countersink Ø6.5

41237 Bone Profiler Ø5.1mm, SP

41238 Bone Profiler Ø6.6mm, SP

41239 Bone Profiler Ø4.1 mm, NP

41240 Bone Remover Cover Screw, SP

41241 Bone Remover Cover Screw, NP

Ø	19mm
Ø	2.3mm
Ø	2.8mm
Ø	3.3mm
Ø	3.8mm
Ø	4.3mm

Ø	2.7mm
Ø	2.9mm
Ø	3.3mm
Ø	3.5mm
Ø	4.3mm
Ø	4.7mm
Ø	2.0mm
Ø	2.8mm
Ø	3.3mm
Ø	3.7mm
Ø	4.2mm
Ø	1.8mm
Ø	2.6mm
Ø	3.0mm
Ø	3.5mm
Ø	4.8mm
Ø	5.2mm

Multiple Use

Kits

Surgical & Prosthetic Trays

Art. No. Description

51180	Neoss Prosthetic Tray – Neo
51186	Neoss ProActive® Implant System Box
51189	Neoss Surgical and Prosthetic Tray – ProActive® Edge
51188	Neoss Surgical and Prosthetic Tray – ProActive® Tapered
51187	Neoss Surgical and Prosthetic Tray – ProActive® Straight



Neoss Prosthetic Tray – Neo



Neoss ProActive® Implant System Box



ProActive® Edge



ProActive® Tapered



ProActive® Straight

Surgical & Prosthetic Kits with Instruments and Drills

Art. No. Description

51181	Neoss Prosthetic Instrument Kit – Neo <i>Includes art no 51121, 51126, 51139, 51140, 51141 and 51180</i>
41151	Surgical Instrument Kit <i>Includes art no 41120, 51118, 51119, 51125, 51126, 51139, 51140 and 51144</i>
41254	Drill and Surgical Instrument Kit, Edge Implants Ø3.5 – 5.0 <i>Single use twist drills. Includes art no 41151, 41152, 41153, 41170, 41177, 41210, 41211, 41212, 41213, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51189 and 51190</i>
41253	Drill and Surgical Instrument Kit, Edge Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41151, 41152, 41153, 41170, 41210, 41211, 41212, 41213, 41236, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51189 and 51190</i>
41252	Drill and Surgical Instrument Kit, Tapered Implants Ø3.5 – 5.0 <i>Single use twist drills. Includes art no 41151, 41152, 41153, 41158, 41160, 41162, 41164, 41170, 41178, 41179, 41180, 41181, 41192, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51188 and 51190</i>
41251	Drill and Surgical Instrument Kit, Tapered Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41151, 41152, 41153, 41170, 41178, 41179, 41180, 41181, 41218, 41220, 41222, 41224, 41235, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51188 and 51190</i>
41248	Drill and Surgical Instrument Kit, Straight Implants Ø3.5 – 5.0 <i>Single use twist drills. Includes art no 41110, 41111, 41112, 41133, 41151, 41152, 41153, 41158, 41160, 41162, 41164, 41170, 41193, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51187 and 51190</i>
41247	Drill and Surgical Instrument Kit, Straight Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41110, 41111, 41112, 41133, 41151, 41152, 41153, 41170, 41218, 41220, 41222, 41224, 41233, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51187 and 51190</i>
41250	Drill and Surgical Instrument Kit, Straight Implants NP Ø3.25 <i>Single use twist drills. Includes art no 41120, 41139, 41152, 41153, 41154, 41156, 41157, 41170, 41239, 41241, 51121, 51125, 51126, 51137, 51139, 51140, 51141, 51142, 51144, 51145, 51186, 51187 and 51190</i>
41249	Drill and Surgical Instrument Kit, Straight Implants NP Ø3.25 <i>Multiple use twist drills. Includes art no 41120, 41139, 41152, 41153, 41170, 41215, 41216, 41217, 41239, 41241, 51121, 51125, 51126, 51137, 51139, 51140, 51141, 51142, 51144, 51145, 51186, 51187 and 51190</i>
41256	Drill and Surgical Instrument Add-On-Kit, Straight Implant Ø3.25 <i>Single use twist drills. Includes art no 41139, 41156, 41239, 41241, 51142 and 51145</i>
41255	Drill and Surgical Instrument Add-On-Kit, Straight Implant Ø3.25 <i>Multiple use twist drills. Includes art no 41139, 41216, 41239, 41241, 51142 and 51145</i>



Neoss Prosthetic Instrument Kit – Neo



Drill and Surgical Instrument Kit, Edge



Drill and Surgical Instrument Kit, Tapered



Drill and Surgical Instrument Kit, Straight



Drill and Surgical Instrument Kits also includes the prosthetic instruments and a Depth Gauge Probe. All in one box for your convenience.

Neoss® Guide Kits Ø3.25 – Ø5.0 Hard Stop & Flexible



Neoss Guide Kit
Ø3.25 – Ø5.0 G5 Hard Stop



Neoss Guide Kit
Ø3.25 – Ø5.0 G3/5/7 Flexible

Art. No.	Description
51200	Neoss Guide Kit Ø3.25 – Ø5.0 G5 Hard Stop
51258	Neoss Guide Kit Ø3.25 – Ø5.0 G3/5/7 Flexible
51260	Neoss Guide Instrument Kit <i>Refer to page 15</i>

Included in both kits (51200 & 51258)

Art. No.	Description
51201	Neoss Guide – Punch – Narrow
51203	Neoss Guide – Start Drill – Narrow
51204	Neoss Guide – Punch – Wide
51206	Neoss Guide – Start Drill – Wide
51207	Neoss Guide – Drill Hub I, Narrow
51208	Neoss Guide – Drill Hub II, Narrow
51209	Neoss Guide – Drill Hub III, Narrow
51210	Neoss Guide – Drill Hub IV, Narrow
51211	Neoss Guide – Drill Hub V, Narrow
51212	Neoss Guide – Drill Hub I, Wide
51213	Neoss Guide – Drill Hub II, Wide
51214	Neoss Guide – Drill Hub III, Wide
51215	Neoss Guide – Drill Hub IV, Wide
51216	Neoss Guide – Drill Hub V, Wide
51217	Neoss Guide – Twist Drill, Ø2.2 × 7 – 15 mm
51218	Neoss Guide – Twist Drill, Ø2.85 × 7 – 15 mm
51219	Neoss Guide – Twist Drill, Ø3.0 × 7 – 15 mm
51220	Neoss Guide – Twist Drill, Ø3.2 × 7 – 15 mm
51221	Neoss Guide – Twist Drill, Ø3.4 × 7 – 15 mm
51222	Neoss Guide – Twist Drill, Ø3.6 × 7 – 15 mm
51223	Neoss Guide – Twist Drill, Ø3.9 × 7 – 15 mm
51224	Neoss Guide – Twist Drill, Ø4.1 × 7 – 15 mm
51225	Neoss Guide – Twist Drill, Ø4.4 × 7 – 15 mm
51226	Neoss Guide – Twist Drill, Ø4.6 × 7 – 15 mm
51238	Neoss Guide – Implant Inserter, SP, Wide
51239	Neoss Guide – Implant Inserter, SP, Narrow
51240	Neoss Guide – Implant Inserter, NP, Narrow
51241	Neoss Guide Kit System Box



Punch, Narrow



Punch, Wide



Start Drill, Narrow





Start Drill, Wide



Drill Hub, Narrow



Drill Hub, Wide



Twist Drill



Implant Inserter, SP, Narrow




Implant Inserter, SP, Wide



Implant Inserter, NP, Narrow

Hard Stop (51200)

Art. No.	Description
51246	Neoss Guide – Mill – Narrow G5
51247	Neoss Guide – Mill – Wide G5
51249	Neoss Guide – Countersink, Ø3.25 S – Narrow G5
51250	Neoss Guide – Countersink, Ø3.5 S/T/E – Narrow G5
51251	Neoss Guide – Countersink, Ø3.5 S/T/E – Wide G5
51252	Neoss Guide – Countersink, Ø4.0 S – Wide G5
51253	Neoss Guide – Countersink, Ø4.5 S – Wide G5
51254	Neoss Guide – Countersink, Ø5.0 S – Wide G5
51255	Neoss Guide – Countersink, Ø4.0 T/E – Wide G5
51256	Neoss Guide – Countersink, Ø4.5 T/E – Wide G5
51257	Neoss Guide – Countersink, Ø5.0 T/E – Wide G5
51261	Neoss Guide add-on-kit, Countersinks and Mill G3/5/7 Flexible (incl. 51202, 51205, 51227, 51228, 51248, 51229, 51230, 51231, 51232, 51233, 51234)



Flexible (51258)


Art. No.	Description
51202	Neoss Guide – Mill – Narrow G3/5/7
51205	Neoss Guide – Mill – Wide G3/5/7
51227	Neoss Guide – Countersink, Ø3.25 S – Narrow G3/5/7
51228	Neoss Guide – Countersink, Ø3.5 S/T/E – Narrow G3/5/7
51248	Neoss Guide – Countersink, Ø3.5 S/T/E – Wide G3/5/7
51229	Neoss Guide – Countersink, Ø4.0 S – Wide G3/5/7
51230	Neoss Guide – Countersink, Ø4.5 S – Wide G3/5/7
51231	Neoss Guide – Countersink, Ø5.0 S – Wide G3/5/7
51232	Neoss Guide – Countersink, Ø4.0 T/E – Wide G3/5/7
51233	Neoss Guide – Countersink, Ø4.5 T/E – Wide G3/5/7
51234	Neoss Guide – Countersink, Ø5.0 T/E – Wide G3/5/7
51259	Neoss Guide add-on-kit, Countersink and Mill G5 Hard Stop (incl. 51246, 51247, 51249, 51250, 51251, 51252, 51253, 51254, 51255, 51256, 51257)



Assortments

51260 includes: Instruments

Art. No.	Description
51121	Ratchet
51144	Wrench Adapter
51137	Implant Insert, SP 17mm
51126	Manual Handle
51139	Neo Screwdriver, Machine 22mm



Guide Sleeves

Guide Sleeves can be ordered separately for Guide.


Art. No.	Description
51242	Neoss Guide – Open Sleeve, Narrow – 5 pcs
51243	Neoss Guide – Closed Sleeve, Narrow – 5 pcs
51244	Neoss Guide – Open Sleeve, Wide – 5 pcs
51245	Neoss Guide – Closed Sleeve, Wide – 5 pcs



Implant Mounts (Optional)

Implant Mounts are optional and can be ordered separately.

Art. No.	Description
51235	Neoss Guide – Implant Mount, NP, Narrow G5
51236	Neoss Guide – Implant Mount, SP, Narrow G5
51237	Neoss Guide – Implant Mount, SP, Wide G5



Neoss Guide add-on kits allow you to add components or convert between workflows, with 51261 providing Countersinks and Mill G3/5/7 for flexible mode, and 51259 providing Countersink and Mill G5 for hard stop mode. Please contact your local Neoss representative for available options.

Instruments

Replacement components for legacy guide kits still available – please contact your local Area Manager.

Art. No. Description

51137	Implant Inserter, SP 17 mm
51118	Implant Inserter, SP 22 mm
51119	Implant Inserter, SP 32 mm
51145	Implant Inserter, NP 24 mm
51142	Implant Inserter, NP 32 mm
51139	Neo Screwdriver, Machine 22 mm
51140	Neo Screwdriver, Machine 32 mm
51141	Neo Screwdriver, Manual 15 mm
51183	iGO Screwdriver, Machine 22 mm
51184	iGO Screwdriver, Machine 32 mm
51125	Direction/Depth Gauge – 4 pcs
41120	Drill Extender
51190	Depth Gauge Probe
51121	Ratchet
51126	Manual Handle
51191	Insertion Handle – Large
51144	Wrench Adapter
51149	Angulation Gauge, SP
90138	Adjustment Handle
41214	X-ray Planner ProActive Edge Implants
41119	X-ray Planner ProActive Straight Implants
41183	X-ray Planner ProActive Tapered Implants
41203	X-ray Planner Ø6.0 Implants
41200	X-ray Planner Ø6.5 Implants
51156	Neoss® Guide – Sleeve Ø5.0 – 10 pcs
51157	Neoss® Guide – Sleeve Ø6.0 – 5 pcs
51158	Drill Sleeve Ø2.2 – 10 pcs

NEO

iGO




Demo Kit

Art. No. Description

61150	King Size Dummy Kit <i>Includes Neoss Straight, Tapered and Edge Implants, Implant Inserter, Screwdriver, Abutment Screw, NeoLink® Abutment, Esthetic Healing Abutment Narrow Incisor with ScanPeg and Crown</i>
-------	---



The *measurement of implant stability* begins with NeoTell®



Resonance Frequency Device (RFA)

NeoTell®

Art. No. Description

91060	NeoTell® Instrument
<i>Includes NeoTell® Instrument, SmartPeg Mount, TestPeg and USB cable, type A – C</i>	
91070	SmartPeg Type 18 – Neoss SP (5 pcs)
91071	SmartPeg Type 38 – Neoss NP (5 pcs)
91072	SmartPeg Type 52 – Neoss Access Abutment (5 pcs)
91073	SmartPeg Type 98 – Neoss Multi Unit Abutment (5 pcs)
91075	Neoss SmartPeg Mount

SmartPegs for other brands are available, contact your local representative for details

New opportunities for your *EasyDigital* workflow



NeoScan™ 1000

Intraoral Scanners

NeoScan™ 1000

Art. No. Description

91000	NeoScan™ 1000 Intraoral Scanner
<i>Includes Handpiece, 4 tips (2×S, 2×L), Holder Set, Manuals, Suitcase. Extended warranties are available for purchase</i>	
91012	NeoScan™ 1000 – Tip L TP102N (2 pcs)
91013	NeoScan™ 1000 – Tip S TP103N (2 pcs)
91022	Trolley for Intraoral Scanner
91024	NeoScan™ 1000 – Warranty extension – 1 year
91026	NeoScan™ 1000 – Warranty extension – 2 years
91027	NeoScan™ 1000 – Warranty extension – 3 years




NeoScan™ 2000

NeoScan™ 2000

Art. No. Description

91001	NeoScan™ 2000 Intraoral Scanner
<i>Includes Handpiece, 4 tips (2×S, 2×L), Charging holder, 2 batteries, Stand, Manuals, Wi-Fi adaptor, Suitcase. Extended warranties are available for purchase</i>	
91014	NeoScan™ 2000 – Tip L TP102N (2 pcs)
91015	NeoScan™ 2000 – Tip S TP103N (2 pcs)
91040	NeoScan™ 2000 – Battery (2 pcs)
91032	NeoScan™ 2000 – Warranty extension – 1 year
91033	NeoScan™ 2000 – Warranty extension – 2 years
91034	NeoScan™ 2000 – Warranty extension – 3 years



Digital Impression

Esthetic Healing Abutments with ScanPeg

All Esthetic Healing Abutments include Titanium Screw and a ScanPeg

Art. No.	Description	A	B	C
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor	5 mm	5.7 mm	6.1 mm
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor		5.2 mm	5 mm
31362	Esthetic Healing Abutment with ScanPeg – Canine		6.3 mm	5.4 mm
31363	Esthetic Healing Abutment with ScanPeg – Pre-molar		6.8 mm	5.8 mm
31364	Esthetic Healing Abutment with ScanPeg – Molar		7.8 mm	7.3 mm

Esthetic Healing Abutment with ScanPeg



Esthetic Healing Abutments

Scan Bodies

Art. No.	Description
31451	Scan Body Ti Neoss® NP – 1 pc
31452	Scan Body Ti Neoss® SP – 1 pc
31453	Scan Body Ti Neoss® Access – 1 pc
31458	Scan Body Ti Neoss® Multi-Unit – 1 pc
31459	Scan Body Ti Neoss® NP – 5 pcs
31460	Scan Body Ti Neoss® SP – 5 pcs
31461	Scan Body Ti Neoss® Access – 5 pcs
31462	Scan Body Ti Neoss® Multi-Unit – 5 pcs



Scan Body Ti Neoss® Access Scan Body Ti Neoss® Multi-Unit

ScanPeg


Art. No.	Description
31365	ScanPeg – 5 pcs Used only with Esthetic Healing Abutments



ScanPeg

ScanKey

Art. No.	Description
	ScanKey Kit The ScanKey Kits are available as STL files for lab printing. Download the ScanKey at neoss.com.
	ScanKey, NP The ScanKey Kit – NP is compatible with the Scan Body 31452, 31453 and 31458.
	ScanKey, Multi-Unit The ScanKey Kit – Multi-Unit is compatible with the Scan Body 31451, 31459 and 31462.




ScanKey Kit– 6 pcs

Digital Analogs

Art. No.	Description
31470	Digital Analog Neoss® NP – 1 pc
31471	Digital Analog Neoss® SP – 1 pc
31472	Digital Analog Neoss® Access – 1 pc
31473	Digital Analog Neoss® Multi-Unit – 1 pc
31474	Digital Analog Neoss® NP – 5 pcs
31475	Digital Analog Neoss® SP – 5 pcs
31476	Digital Analog Neoss® Access – 5 pcs
31477	Digital Analog Neoss® Multi-Unit – 5 pcs
31478	Digital Analog Neoss® Fixation Screw – 1 pc
31479	Digital Analog Neoss® Fixation Screw – 5 pcs
	Digital Analog Neoss® Fixation Tool The Fixation Tool is available as an STL file for lab printing. Download the Tool STL file at neoss.com.





Digital Prosthetics

NeoBase® Abutments SSC (Straight Screw Channel)*

also referred to as TiBase

All NeoBase® Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw



NeoBase® Abutments SSC



NeoBase® Abutments SSC, NP

Art. No.	Description	A	B	C
31393	NeoBase® Mono SSC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31394	NeoBase® Mono SSC G1.5mm – H3.6mm	5.1mm	1.5mm	
31395	NeoBase® Mono SSC G0.3mm – H5.6mm	5.9mm	0.3mm	
31396	NeoBase® Mono SSC G1.5mm – H5.6mm	7.1mm	1.5mm	
31397	NeoBase® Multi SSC G0.3mm – H3.6mm	3.9mm	0.3mm	3.8mm
31398	NeoBase® Mono SSC G0.3mm – H3.6mm, NP	3.9mm		
31399	NeoBase® Mono SSC G0.3mm – H5.6mm, NP	5.9mm	0.3mm	
31400	NeoBase® Multi SSC G0.3mm – H3.6mm, NP	3.9mm		

Includes Neo Abutment Screw NP and Neo Laboratory Screw NP

*See page 26 for ARC offering – custom-made implant prosthetics

NeoBase® Abutments ASC (Angulated Screw Channel)*

also referred to as TiBase

All NeoBase® Abutments ASC include iGO Abutment Screw and iGO Laboratory Screw



NeoBase® Abutments ASC



NeoBase® Abutments ASC, NP

Art. No.	Description	A	B	C
31367	NeoBase® Mono ASC G 0.3 mm – H 3.6 mm	3.9 mm	0.3 mm	4.5 mm
31368	NeoBase® Mono ASC G 1.5 mm – H 3.6 mm	5.1 mm	1.5 mm	
31369	NeoBase® Mono ASC G 0.3 mm – H 5.6 mm	5.9 mm	0.3 mm	
31370	NeoBase® Mono ASC G 1.5 mm – H 5.6 mm	7.1 mm	1.5 mm	
31371	NeoBase® Multi ASC G 0.3 mm – H 3.6 mm	3.9 mm	0.3 mm	3.8 mm
31372	NeoBase® Mono ASC G 0.3 mm – H 3.6 mm, NP	3.9 mm	0.3 mm	
31373	NeoBase® Mono ASC G 0.3 mm – H 5.6 mm, NP	5.9 mm		
31374	NeoBase® Multi ASC G 0.3 mm – H 3.6 mm, NP	3.9 mm		

Includes iGO Abutment Screw NP and iGO Laboratory Screw NP

*See page 26 for ARC offering – custom-made implant prosthetics

Neoss TiBase Abutments (compatible with Dentsply Sirona®)


All TiBase Abutments include Neo Abutment Screw and Titanium Screw



TiBase Abutments

Neoss ScanPost

Art. No.	Description	A	B	C
31329	Neoss TiBase N (NB B 3.4 L)	5.3mm	0.6mm	4.5mm
31330	Neoss TiBase W (NB B 4.1 L)			5.0mm
31345	Neoss TiBase, NP (FX 3.4 S)	5mm	0.3mm	3.9mm
31331	Neoss ScanPost L (TiBase) Includes ScanPost with Screw			
31346	Neoss ScanPost S (TiBase, NP) Includes ScanPost with Screw NP			



Ti NeoLink®

Titanium NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings



Ti NeoLink®

Art. No.	Description	A	B	C
31133	Ti NeoLink® Mono	3.9mm	0.3mm	4.5mm
31134	Ti NeoLink® Multi			3.8mm
31234	Ti NeoLink® Mono, NP			

Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings

Impressions & Replicas

Impression Copings


Art. No.	Description
31170	Impression Coping 8mm and Replica
31145	Impression Coping 11mm and Replica
31183	Impression Coping 18mm and Replica <i>31170, 31145 and 31183 include Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31302	Impression Coping Open Tray 8mm and Replica
31185	Impression Coping Open Tray 13mm and Replica <i>31302 and 31185 include Coping, Screw, Replica and Extension Tube</i>
31197	Impression Coping 11mm and Replica, NP <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31442	Multi-Unit Impression Coping Open tray 6mm*
31443	Multi-Unit Impression Coping Closed tray 10mm*
31366	Access Impression Coping 13mm and Replica
31186	Impression Coping Cap Closed Tray – 10 pcs
31262	Impression Coping Extension Tube – 20 pcs

*Not co-packed and should be ordered separately


Replicas

Art. No.	Description
31444	Multi-Unit Replica – 5 pcs*
31231	Access Replica – 5 pcs
31146	Implant Replica – 5 pcs
31198	Implant Replica, NP – 5 pcs
31304	Implant Replica, NP – 1 pc
31166	Protection Replica – 1 pc


*Not co-packed and should be ordered separately



Impression Coping and Replica



Impression Coping Open Tray and Replica



Impression Coping and Replica, NP



Access Impression Coping and Replica



Impression Coping Cap

Impression Coping Extension Tube



Multi-Unit Replica



Access Replica




Implant Replica, SP



Implant Replica, NP



Protection Replica



Final Abutments

Preparable Titanium Abutments

All Preparable Abutments include Neo Abutment Screw and Neo Laboratory Screw



Art. No.	Description	A	B	C	D	E
31280	Preparable Ti Abutment 1.5mm 0° – Wide Incisor	8.5mm	2.5mm	1.5mm	5.8mm	6.5mm
31281	Preparable Ti Abutment 1.5mm 15° – Wide Incisor					
31282	Preparable Ti Abutment 1mm 0° – Narrow Incisor	8mm	2mm	1mm	5mm	5.1mm
31283	Preparable Ti Abutment 1mm 15° – Narrow Incisor					
31284	Preparable Ti Abutment 3mm 0° – Narrow Incisor	10mm	4mm	3mm	5.1mm	5.1mm
31285	Preparable Ti Abutment 3mm 15° – Narrow Incisor					
31286	Preparable Ti Abutment 3mm 0° – Canine	10mm	4mm	3mm	6.7mm	5.4mm
31305	Preparable Ti Abutment 1.5mm 0° – Pre-molar	7mm	2mm	1.5mm	6.9mm	5.7mm
31287	Preparable Ti Abutment 3mm 0° – Pre-molar	9mm	3mm	3mm	6.9mm	5.7mm
31306	Preparable Ti Abutment 1.5mm 0° – Molar	7mm	2mm	1.5mm	7.8mm	7.1mm
31288	Preparable Ti Abutment 3mm 0° – Molar	9mm	3mm	3mm	7.8mm	7.1mm
31201	Preparable Ti Abutment 1mm 0°, NP	8mm	2mm	1mm	4mm	
31202	Preparable Ti Abutment 1mm 15°, NP					

Preparable Titanium Abutments, Alternative Emergence Profiles




Art. No.	Description	A	B	C	D	E
31317	Preparable Ti Abutment 2mm 0° – Wide Emergence <i>Includes Neo Abutment Screw, Neo Laboratory Screw, specific Cover Screw and specific Healing Abutment PEEK with Screw. Compatible with Implant diameters Ø5.0 – 6.5mm</i>	9mm	2mm	2mm	7.9mm	7.1mm
31318	Preparable Ti Abutment 1.5mm 0° – Narrow Emergence	8.5mm	2mm	1.5mm	4.4mm	4.5mm
31319	Preparable Ti Abutment 1.5mm 15° – Narrow Emergence <i>All Narrow Emergence Preparable Abutments include Neo Abutment Screw and Neo Laboratory Screw. Only use with Ø3.5 and Ø4.0 mm implants</i>					

Zirconia Abutments

All Zirconia Abutments include separate Zirconia Coping, Ti NeoLink® Mono, Neo Abutment Screw and Neo Laboratory Screw



Art. No.	Description	A	B	C	D	E
31290	Zirconia Abutment – Wide Incisor				5.8mm	6.4mm
31291	Zirconia Abutment – Narrow Incisor		2.5mm	1.5mm	5.4mm	5.1mm
31292	Zirconia Abutment – Canine	9mm			6.2mm	5.8mm
31293	Zirconia Abutment – Pre-molar		2mm	2mm	6.4mm	5.5mm
31294	Zirconia Abutment – Molar				7.2mm	6.7mm




Final Restorations – Casting

Gold NeoLink® (Direct Casting)

Gold NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings


Art. No.	Description	A	B	C
31102	Gold NeoLink® Mono	4.0 mm	0.4 mm	4.5 mm
31103	Gold NeoLink® Multi			3.8 mm
31199	Gold NeoLink® Mono, NP <i>Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings</i>			



CoCr Abutments (Indirect Casting)

CoCr Abutments include Neo Abutment Screw, Neo Laboratory Screw and Straight Coping


Art. No.	Description	A	B	C
31269	CoCr Abutment Mono	15.7 mm	0.7 mm	4.5 mm
31260	CoCr Abutment Multi			



Burnout Abutments (Direct Casting)

31333, 31157 and 31266 include Straight Coping

Art. No.	Description	A	B	C
31333	Burnout Abutment Mono – 2 pcs	15.7 mm	0.7 mm	4.5 mm
31157	Burnout Abutment Multi – 2 pcs			3.8 mm
31266	Access Burnout Abutment Multi – 2 pcs			
31255	Burnout NeoLink® Multi NP – 2 pcs <i>Includes set of Preformed Plastic Copings</i>	4.3 mm		







Neoss® Medentika Blanks

Each 5 – pcs includes 5 laboratory and 5 Prosthetics Screws
Available in Titanium and Cobalt Chrome

Art. No.	Description
31488	Neoss® Medentika Blank Ti SSC, Ø11.5 mm, SP – 5 pcs
31489	Neoss® Medentika Large Blank Ti SSC, Ø16 mm, SP – 5 pcs
31490	Neoss® Medentika Blank Ti SSC, Ø11.5 mm, NP – 5 pcs
31491	Neoss® Medentika Blank CoCr SSC, Ø11.5 mm, SP – 5 pcs
31492	Neoss® Medentika Blank CoCr SSC, Ø11.5 mm, NP – 5 pcs



Experience *easy-to-use full arch flexibility* with predictable results using Neoss4+ treatment solution



Neoss4+ treatment solution

Multi-Unit Abutments


Clinical screw fixed within the Multi-Unit Abutment
Healing Cap and Replica need to be ordered separately

Art. No.	Description	A	B	C
31411	Multi-Unit Abutment 1.0mm	1mm	4.6mm	
31412	Multi-Unit Abutment 2.0mm	2mm		
31413	Multi-Unit Abutment 3.0mm	3mm		
31414	Multi-Unit Abutment 4.0mm	4mm		
31425	Multi-Unit Abutment 10° 2.0mm	2mm	4.6mm	
31426	Multi-Unit Abutment 10° 3.0mm	3mm		
31427	Multi-Unit Abutment 10° 4.0mm	4mm		
31429	Multi-Unit Abutment 17° 2.5mm	2.5mm	4.6mm	
31430	Multi-Unit Abutment 17° 3.5mm	3.5mm		
31431	Multi-Unit Abutment 17° 4.5mm	4.5mm		
31433	Multi-Unit Abutment 30° 3.0mm	3mm	4.6mm	
31434	Multi-Unit Abutment 30° 4.0mm	4mm		
31435	Multi-Unit Abutment 30° 5.0mm	5mm		
31440	Multi-Unit Healing Cap	5mm	5mm	
31442	Multi-Unit Impression Coping Open tray 6mm*			
31443	Multi-Unit Impression Coping Closed tray 10mm*			
31444	Multi-Unit Replica – 5 pcs			
31445	Multi-Unit Provisional Abutment Ti Includes Provisional/Laboratory Screw	15.7mm	0.7mm	5mm
31446	Multi-Unit NeoBase® SSC G 0.5mm – H 4.5mm	5mm	0.5mm	5.1mm
31447	Multi-Unit NeoBase® ASC G 0.5mm – H 4.5mm iGO	5mm	0.5mm	5.1mm
31458	Scan Body Ti Neoss® Multi-Unit – 1 pc	See page 18 for details		
31462	Scan Body Ti Neoss® Multi-Unit – 5 pcs			
31473	Digital Analog Neoss® Multi-Unit – 1 pc	See page 18 for details		
31477	Digital Analog Neoss® Multi-Unit – 5 pcs			


*Not co-packed and should be ordered separately

Multi-Unit/Access Screws

Art. No.	Description
31229	Multi-Unit/Access Neo Prosthetic Screw – 1pc
31257	Multi-Unit /Access Neo Provisional/Laboratory Screw – 5 pcs
31265	Multi-Unit/Access Neo Laboratory Screw Long – 5 pcs
31348	Multi-Unit/Access iGO Prosthetic Screw – 1pc
31351	Multi-Unit/Access iGO Laboratory Screw – 1pc



A=Height
B=Diameter



A=Height
B=Margin height
C=Margin diameter

Access Abutments

Clinical screw fixed within the Multi-Unit Abutment

Art. No. Description

A B

All Access Abutments include Abutment and a Healing Abutment

31238	Access Abutment 1.5mm	1.5mm	4mm
31215	Access Abutment 3.0mm	3mm	
31216	Access Abutment 4.0mm	4mm	
31218	Access Abutment 6.0mm	6mm	
31219	Access Abutment 10° 2.6mm	2.6mm	4mm
31243	Access Abutment 10° 4.6mm	4.6mm	
31221	Access Abutment 20° 2.6mm	2.6mm	
31244	Access Abutment 20° 4.6mm	4.6mm	
31223	Access Abutment 30° 2.9mm	2.9mm	
31245	Access Abutment 30° 4.9mm	4.9mm	
31225	Access Impression Coping 8mm and Replica <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>		
31366	Access Impression Coping 13mm and Replica		
31453	Scan Body Ti Neoss® Access - 1 pc	See page 18 for details	
31461	Scan Body Ti Neoss® Access - 5 pcs	See page 18 for details	
31472	Digital Analog Neoss® Access - 1 pc	See page 18 for details	
31476	Digital Analog Neoss® Access - 5 pcs	See page 18 for details	



Access Abutments



Access Impression Coping and Replica



Scan Body Ti Neoss® Access



Digital Analog Neoss® Access




Access Ti NeoLink®



Access Gold NeoLink®



Access Equator Abutments



A=Height
B=Margin height
C=Margin diameter

Art. No. Description


A B C

31226	Access Ti NeoLink® Multi	4.3mm	0.3mm	4.5mm
31227	Access Gold NeoLink® Multi <i>31226 and 31227 include NeoLink®, Prosthetic Screw, Provisional/Laboratory Screw and a Plastic Coping</i>	4.7mm	0.7mm	4.5mm
90207	Access Equator Abutment Kit 0.5mm			
90208	Access Equator Abutment Kit 1.5mm			
31301	Access Provisional Ti Abutment Multi <i>Includes Provisional/Laboratory Screw and Straight Coping</i>	15.7mm	0.7mm	4.5mm
31266	Access Burnout Abutment Multi – 2pcs <i>Includes Straight Coping</i>			

Multi-Unit/Access Wide Screws

Art. No. Description

31480	Neo Provisional/Laboratory Wide Screw, Multi-Unit/Access - 1 pc
31481	Neo Provisional/Laboratory Wide Screw, Multi-Unit/Access - 5 pcs
31482	Neo Prosthetic Wide Screw, Multi-Unit/Access - 1 pc
31483	Neo Prosthetic Wide Screw, Multi-Unit/Access - 5 pcs
31484	iGO Provisional/Laboratory Wide Screw, Multi-Unit/Access - 1 pc
31485	iGO Provisional/Laboratory Wide Screw, Multi-Unit/Access - 5 pcs
31486	iGO Prosthetic Wide Screw, Multi-Unit/Access - 1 pc
31487	iGO Prosthetic Wide Screw, Multi-Unit/Access - 5 pcs



Access Provisional Ti Abutment Multi

Access Burnout Abutment Multi


Coping



iGO Provisional/Laboratory Wide Screw, Multi-Unit/Access



iGO Prosthetic Wide Screw, Multi-Unit/Access



Neo Provisional/Laboratory Wide Screw, Multi-Unit/Access



Neo Prosthetic Wide Screw, Multi-Unit/Access


Overdenture Abutments

Equator Abutments

All kits include Equator Abutment, Housing with processing Cap, Protector Disk, Standard, Strong, Soft and Extra Soft Retention caps

Art. No. Description

90201	Neoss Equator Abutment Kit 1mm	90210	Neoss Equator Cap and Housing Kit
90202	Neoss Equator Abutment Kit 2mm	90215	Neoss Equator Cap and Housing Kit – Extended Divergence
90242	Neoss Equator Abutment Kit 2mm, NP	90211	Neoss Equator Cap Standard Retention, White – 2 pcs
90203	Neoss Equator Abutment Kit 3mm	90212	Neoss Equator Cap Strong Retention, Violet – 2 pcs
90243	Neoss Equator Abutment Kit 3mm, NP	90213	Neoss Equator Cap Soft Retention, Pink – 2 pcs
90204	Neoss Equator Abutment Kit 4mm	90214	Neoss Equator Cap Extra Soft Retention, Yellow – 2 pcs
90244	Neoss Equator Abutment Kit 4mm, NP	90220	Neoss Equator Impression Coping & Replica – 2 sets
90205	Neoss Equator Abutment Kit 5mm	90222	Neoss Equator Protector Disk – 10 pcs
90206	Neoss Equator Abutment Kit 6mm	90232	Neoss Equator Cap Tool
90207	Access Equator Abutment Kit 0.5mm		
90208	Access Equator Abutment Kit 1.5mm		



Screw Guide

Screws

Art. No. Description

31258	Neo Abutment Screw, SP – 1pc	31484	iGO Provisional/Laboratory Wide Screw, Multi-Unit – 1pc
31259	Neo Abutment Screw, SP – 5pcs	31485	iGO Provisional/Laboratory Wide Screw, Multi-Unit – 5pcs
31240	Neo Abutment Screw Ti, NP – 1pc	31486	iGO Prosthetic Wide Screw, Multi-Unit – 1pc
31149	Neo Provisional/Laboratory Screw, SP – 5pcs	31487	iGO Prosthetic Wide Screw, Multi-Unit – 5pcs
31239	Neo Provisional/Laboratory Screw, NP – 1pc	31229	Multi-Unit/Access Neo Prosthetic Screw – 1pc
31303	Neo Laboratory Screw Long, SP – 5pcs	31257	Multi-Unit/Access Neo Provisional/Laboratory Screw – 5pcs
31347	iGO Abutment Screw, SP – 1pc	31265	Multi-Unit/Access Neo Laboratory Screw Long – 5pcs
31349	iGO Abutment Screw, NP – 1pc	31480	Neo Provisional/Laboratory Wide Screw, Multi-Unit – 1pc
31348	Multi-Unit/Access iGO Prosthetic Screw – 1pc	31481	Neo Provisional/Laboratory Wide Screw, Multi-Unit – 5pcs
31350	iGO Laboratory Screw, SP – 1pc	31482	Neo Prosthetic Wide Screw, Multi-Unit – 1pc
31352	iGO Laboratory Screw, NP – 1pc	31483	Neo Prosthetic Wide Screw, Multi-Unit – 5pcs
31351	Multi-Unit/Access iGO Laboratory Screw – 1pc		

Overview abutment, prosthetic and laboratory screws


	Laboratory and Provisional prosthetics				Definitive prosthetics			
	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges
Straight screw channel 	 31149 Neo Laboratory Screw, SP	 31239 Neo Laboratory Screw, NP	 31257 Multi-Unit/Access Neo Provisional/Laboratory Screw	 31480 Neo Provisional/Laboratory Wide Screw, Multi-Unit/Access	 31258 Neo Abutment Screw, SP	 31240 Neo Abutment Screw Ti, NP	 31229 Multi-Unit/Access Neo Prosthetic Screw	 31482 Neo Prosthetic Wide Screw, Multi-Unit/Access
Angulated screw channel 	 31350 iGO Laboratory Screw, SP	 31352 iGO Laboratory Screw, NP	 31351 Multi-Unit/Access iGO Provisional/Laboratory Screw	 31484 iGO Provisional/Laboratory Wide Screw, Multi-Unit/Access	 31347 iGO Abutment Screw, SP	 31349 iGO Abutment Screw, NP	 31348 Multi-Unit/Access iGO Prosthetic Screw	 31486 iGO Prosthetic Wide Screw, Multi-Unit/Access

ARC Solutions Individualised Prosthetics


ARC Solutions are specialists in the manufacture of custom-made implant prosthetics. ARC Solutions places great value on providing customers with maximum service and delivering reliable products.

Get started and how to order

- 1. Go to www.neoss.com/en-au/arc
- 2. Please ensure that you are on the correct country
- 3. Click ORDER NOW
- 4. Log in or create an account
- 5. Select the required product
- 6. Input order details and upload scan files
- 7. Add the order to the cart
- 8. Proceed to checkout and process the order




Scenario 1



I have scan files and would like ARC Solutions to design and manufacture my prosthetic.

- 1. Create a zip file of your files
- 2. Logon to NEOSS CAD/CAM
- 3. Create order basket
- 4. Upload scan files to order
- 5. Check and confirm design proposal
- 6. Receive custom prosthetic

Scenario 2



I have completed design files and would like ARC Solutions to manufacture my prosthetic.

- 1. Create a zip file of your files
- 2. Logon to NEOSS CAD/CAM
- 3. Create order basket
- 4. Upload design files to order
- 5. Receive custom prosthetic

	Abutment						Bridge				Bar			
	Ti		CoCr		ZrO		Ti		CoCr		Ti		CoCr	
	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC
Neoss SP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Neoss NP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Neoss Access							✓	✓	✓	✓	✓	✓	✓	✓
Neoss Multi-Unit							✓	✓	✓	✓	✓	✓	✓	✓

SSC= Straight Screw Channel
ASC= Angulated Screw Channel

ARC Solutions

ARC Solutions offers an attractive bundle of custom-milled titanium abutments or bars in combination with custom-milled zirconia.



Individual abutments

Customised individual abutments are available in titanium and cobalt chromium, with options for straight or angled screw channels to support optimal functional and aesthetic outcomes.



Implant bridges and bars

Precision-milled implant bridges and bars are produced with straight or angled screw channels, ensuring accurate fit, stability and predictable long-term performance across a wide range of indications.



OneToOne

OneToOne is a custom-designed framework available at implant or abutment level, produced in titanium with straight or angled screw channels. An optional zirconia secondary structure is also available, supporting a streamlined and efficient restorative workflow.



FreeD Abutment

The FreeD Abutment is a revolutionary 3D-printed custom abutment solution, available in titanium or cobalt chromium with straight or angled screw channels. By combining a 3D-printed Ti-Base and abutment, the FreeD system delivers exceptional design freedom, accuracy and consistency.



Other products available from ARC



Partial Titanium Bridges
SSC/ASC



Full Arch Bridges
SSC/ASC



Custom Bars

Scan the QR code
or visit neoss.com/arc
to place your order.



Solutions for your guided bone *regeneration* procedures



NeoGenix XP™ Curved Syringe

Bone Graft Materials

NeoGenix XP™

Art. No.	Description	Volume
64097	NeoGenix XP™ 0.25cc (0.25–1mm particle size) Curved Syringe	0.25 cc
64098	NeoGenix XP™ 0.5cc (0.25–1mm particle size) Curved Syringe	0.50 cc



PTFE Suture

Sutures

Art. No.	Description
32.Z4401.00	Neoss PTFE Suture 60 cm 4/0 – Reverse Cutting 16 mm 3/8 (12pcs)
32.Z4408.00	Neoss PTFE Suture 60 cm 6/0 – Reverse Cutting 12 mm 3/8 (12pcs)




Collagen Flex

Resorbable Membranes

Collagen Membranes

Art. No.	Description	A	B
64050	NeoGen® Collagen Membrane Flex S	20 mm	15 mm
64051	NeoGen® Collagen Membrane Flex M	30 mm	20 mm
64052	NeoGen® Collagen Membrane Flex L <i>NeoGen® Collagen Flex resorbs in 3–4 months</i>	40 mm	30 mm






Non-Resorbable Membranes

Ti-Reinforced PTFE Membranes

Art. No.	Description	A	B	C	D
64010	NeoGen® PTFE Membrane, Ti-Reinforced – S I (Small Interproximal)	29 mm	14 mm	19 mm	4 mm
64011	NeoGen® PTFE Membrane, Ti-Reinforced – M I (Medium Interproximal)	30 mm	19 mm	19 mm	9 mm
64012	NeoGen® PTFE Membrane, Ti-Reinforced – L I (Large Interproximal)	36 mm	21 mm	26 mm	11 mm
64013	NeoGen® PTFE Membrane, Ti-Reinforced – M (Medium)	32 mm	22 mm	23 mm	12 mm
64014	NeoGen® PTFE Membrane, Ti-Reinforced – L (Large)	34 mm	25 mm	24 mm	15 mm
64015	NeoGen® PTFE Membrane, Ti-Reinforced – XL (Extra Large)	46 mm	32 mm	24 mm	15 mm




Ti-Reinforced Membranes

NeoGen® Cape

All NeoGen® Cape Spacer include a NeoGen® Cape Screw

Art. No.	Description	A	B	C	D
64040	NeoGen® Cape PTFE Membrane – M	28 mm	30 mm	16 mm	9 mm
64041	NeoGen® Cape PTFE Membrane – L	28 mm	30 mm	22 mm	9 mm
64080	NeoGen® Cape Spacer 0.7 mm, Neoss SP				
64081	NeoGen® Cape Spacer 1.5 mm, Neoss SP				
64082	NeoGen® Cape Spacer 0.7 mm, Neoss NP				
64083	NeoGen® Cape Spacer 1.5 mm, Neoss NP				



Non-Reinforced PTFE Membranes


Art. No.	Description	A	B
64019	NeoGen® PTFE Membrane, Non-Reinforced	18 mm	12 mm
64020	NeoGen® PTFE Membrane, Non-Reinforced	29 mm	14 mm
64021	NeoGen® PTFE Membrane, Non-Reinforced	34 mm	25 mm




Non-Reinforced Membranes

Fixation Kit


Art. No.	Description
64031	Tack Positioning Instrument
64033	Tack Mallet
64034	Tack and Screw Cassette
64036	Membrane Tack 3 mm – 10 pcs
64037	Membrane Screw 3 mm – 10 pcs












ProActive Edge Implants



ProActive Tapered Implants



ProActive Straight Implants

	ProActive Straight Ø3.25 mm	ProActive Straight Ø3.5 mm	ProActive Straight Ø4.0 mm	ProActive Straight Ø4.5 mm	ProActive Straight Ø5.0 mm	ProActive Straight Ø5.5 mm	Drill Stop
 Ø2.2 S							
 Ø2.85 S							
 Ø3.0 S							
 Ø3.2 S							
 Ø3.4 S							
 Ø3.6 S							
 Ø3.8 S							
 Ø4.1 S							
 Ø4.4 S							
 Ø4.6 S							
 Ø4.9 S							
 Ø5.1 S							
 Countersink Ø4.0 S	Ø3.25 Optional use	Ø3.5 Optional use	Ø4.0 Optional use	Ø4.5 Optional use	Ø5.0 Optional use	Ø5.5 Optional use	
 Screw Tap Ø4.0 S	Ø3.25 Optional use	Ø3.5 Optional use	Ø4.0 Optional use	Ø4.5 Optional use	Ø5.0 Optional use	Ø5.5 Optional use	

Bone quality



Regular



Dense

Drill Depth Guides


This guide shows an 11 mm implant in relation to a twist drill and depth guide. Please note actual distance to drill tip is 0.8 mm longer than the reference line.

Note: Depth markings on Lance Drill at 3, 5 and 7 mm, and at 7 and 9 mm on Pilot Drill.

Note: The 3 mm tip cannot be used for depth purposes in conjunction with the Twist Drill, Tapered Ø3.0.



ProActive Wide and Sinus Implants

	ProActive Ø6.0 mm	ProActive Ø6.5 mm
 Ø2.2 S	○	○
 Ø3.0 T		
 Ø3.2 S		
 Ø3.4 T	○	○
 Ø3.6 S		
 Ø3.9 T		
 Ø4.1 S		
 Ø4.4 T		
 Ø4.6 S		
 Ø4.9 T	○	○
 Ø5.1 S		○
 Ø5.5 T	○	
 Ø6.0		○
 Countersink	○	○
	Optional use	Optional use

Bone quality


○ Soft
IV & III

○ Regular
II

Clinical Organisers – Surgical Set-Up


ProActive Edge Implants

Surgical and prosthetic setup




ProActive Tapered Implants

Surgical and prosthetic setup



ProActive Straight Implants

Surgical and prosthetic setup



Prosthetic Solution Standard and Narrow Platform, overview

Overview Prosthetic Solutions, SP and NP (Ø3.25 – 6.5 Implants)



		Units		Retention		Implant platform		Level			Material					Angulated screw channel	Esthetiline
		Mono	Multi	Cement retained	Screw retained	SP	NP***	Implant level	Access level	Multi-Unit Abutment level	Ti	ZrO	Gold	CoCr	Polymer		
Provisional Restorations	Provisional Ti Abutment*	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Esthetic Tissue Former	✓		✓	✓	✓		✓							✓		✓
Final Restorations	Prefabricated	Prepable Ti Abutment	✓		✓	✓	✓	✓			✓						✓
		Zirconia Abutment	✓		✓	✓		✓				✓					✓
		Access Abutment (0 – 30°)		✓		✓		✓			✓						
		Multi-Unit Abutment (0 – 30°)		✓		✓		✓			✓						
	Digital	NeoBase® (SSC, ASC**)	✓	✓	✓	✓	✓	✓		✓	✓					✓	✓
		TiBase***	✓		✓	✓	✓	✓			✓						
		Ti NeoLink®	✓	✓	✓	✓	✓	✓	✓		✓						✓
		Blanks	✓	✓	✓	✓	✓	✓			✓			✓			✓
		Individual Prosthetics** central milling	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓		✓	✓
	Casting	Gold NeoLink®	✓	✓		✓	✓	✓	✓				✓				✓
		Burnout Abutment	✓	✓	✓	✓	✓	✓	✓								
		CoCr Abutment	✓	✓	✓	✓		✓						✓			
Overdenture		Equator Abutment		✓			✓	✓	✓	✓							

* Digital Libraries are available for Provisional Ti Abutment except for NP platform.

** Utilise iGO Abutment Screw and Screwdriver for angulated screw channels, optional for straight screw channels on Individual Prosthetics milled by Neoss

***Compatible with Dentsply Sirona®

Torque and Speed recommendation, overview

Neoss Implant System Torque Recommendation (Ncm)

Implants	Screw Taps	Healing			Provisional		Impression	Final restoration					
		Cover Screw	Healing Abutment PEEK and Esthetic Healing Abutment	Healing Abutment Ti	Esthetic Tissue Former	Provisional Ti Abutment	Impression Coping Screw	Neo Abutment Screw	Multi-Unit and Access Abutment Straight & Angulated	Multi-Unit and Neo Access Prosthetic Screw	Ball Abutment including Access level	Equator Abutment	
												Implant level	Access level
45 Max	40 – 45	10 Max	10 Max	10 Max	20	20	10 Max	32*	32	20	20	32	20

* 35 Ncm optional for high-load cases

Final restoration with iGO Screws		
iGO Abutment Screw, SP	iGO Abutment Screw, NP	iGO Prosthetic Screw, Access
32	20	20

Neoss Implant System Drilling/Insertion Speed Recommendation (rpm)

Drills	Screw Taps	Countersinks	Implants	Bone Mill
800 – 2000**	20	800	20 Max	40 Max

** lower speed for larger drills

Symbols, names and abbreviations



SP = Standard Platform Ø3.5 – 6.5 Implants

NP = Narrow Platform Ø3.25 Implants

iGO = utilises iGO Screw and Screwdriver

Neo = Neoss standard Screw/Screwdriver

SSC = Straight Screw Channel

ASC = Angulated Screw Channel

You can now buy your products online!

Products can also be ordered online at **neoss.com**.

Some products may not be regulatory or commercially available for all markets and available in limited quantities. Contact your local representative for availability and prices.

Disclaimer of Liability

Neoss products may only be used according to the manufacturers' instructions and recommendations.

The user of Neoss products should determine their suitability for particular patients and indications.

Neoss Limited disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or practice in the use or placement of the Neoss products.

Product images are not to scale and are subject to change.

All dimensions are for identification purposes only.

Copyrights, design rights and trademarks:

Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorisation of Neoss Limited.

Neoss, 'N' logo, NeoGen, NeoLink, NeoLoc, NeoBase and Neoss ProActive are registered trademarks of Neoss Group. NeoScan is a trademark of Neoss Group.

Elos Accurate is a registered trademark of Elos Medtech AB. LOCATOR® is a registered trademark of ZEST IP Holdings, LLC U.S.



Do you want to *lower your inventory costs* and minimise confusion?

NeossONE® is a solution unique to the Neoss implant system – one prosthetic platform, across three implant designs. Simply put, the same prosthetic components fit every implant.

The benefits are clear: reducing patient treatment time, optimising inventory control and minimised product confusion.

- ✓ Reduce steps, saves time
- ✓ Less inventory and maximum flexibility
- ✓ Minimise confusion and mistakes
- ✓ Cost-effective solution for all



neossONE®



Order your products
24/7 in our Online Stores



Australia
at shop-au.neoss.com



New Zealand
at shop-nz.neoss.com

AUSTRALIA & NEW ZEALAND

Neoss Australia Pty Ltd

G04 Ground Floor North Tower, John Oxley Centre
339 Coronation Drive, Milton Queensland 4064 Australia
T: +61 7 3216 0165
E: info.au@neoss.com

For our Terms & Conditions consult your local website on
www.neoss.com

Contact your local representative for product availability details.
All products are not available in every market.

10622_26 AUS_NZ © Neoss Group, 2026-08. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorisation of Neoss Group. Product images are not necessarily to scale.

neoss.com